

## CRK Monoclonal Antibody

Catalog No: #27127



Package Size: #27127-1 50ul #27127-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	CRK Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3H7-E5-H8
Isotype	IgG2b
Purification	Affinity purified
Applications	E WB
Species Reactivity	Hu D Rb Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human CRK protein fragments expressed in E.coli.
Target Name	CRK
Other Names	p38, CRKII
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

## Application Details

Predicted MW: 34kd

Western blotting: 1:1000

## Images



Western blot detection of CRK in LNCAP,A431

## Background

The Crk-I and Crk-II forms differ in their biological activities. Crk-II has less transforming activity than Crk-I. Crk-II mediates attachment-induced MAPK8 activation, membrane ruffling and cell motility in a Rac-dependent manner.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.