

Annexin V Polyclonal Antibody

Catalog No: #27222

Package Size: #27222-1 50ul #27222-2 100ul #27222-4 25ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Annexin V Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt Mk
Specificity	This antibody detects endogenous levels of ANXA5 and does not cross-react with related proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human ANXA5 protein fragments expressed in E.coli.
Target Name	Annexin V
Other Names	Anchorin CII; Annexin 5; Annexin A5; Annexin V; Annexin-5; Annexin5; AnnexinA5; AnnexinV; ANX 5; ANX A5; ANX5; ANXA5; ANXA5_HUMAN; Calphobindin I; CBP I; CBP-I; Endonexin II; ENX 2; ENX2; Lipocortin V; PAP I; PAP-I; Placental anticoagulant protein 4;
Accession No.	Uniprot: P08758 Gene ID: 308
SDS-PAGE MW	30kd
Formulation	Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Storage	store at -20 Λ C

Application Details

Western blotting: 1:2000

Images



Western blot detection of ANXA5 antibody in 3T3, HeLa, C6, COS7 and MCF7 cell lysates using ANXA5 antibody (1:2000 diluted). Predicted band size: 30 kDa. Observed band size: 30 kDa. Exposure time: 20s.

Background

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.