

JUN Antibody

Catalog No: #32039

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	JUN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total JUN protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human JUN.
Target Name	JUN
Other Names	JUN; AP-1; AP1; c-Jun;
Accession No.	Swiss-Prot:P05412NCBI Gene ID:3725
SDS-PAGE MW	36KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:500 - 1:2000:

Immunofluorescence: 1:50 - 1:200

Images



Western blot analysis of extracts of PC3 cells, using c-Jun antibody.



Immunofluorescence analysis of A549 cell using Jun antibody.
Blue: DAPI for nuclear staining.

Background

c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1). AP-1 is composed of dimers of Fos, Jun and ATF family members and binds to and activates transcription at TRE/AP-1 elements (reviewed in 1). Extracellular signals including growth factors, chemokines and stress activate AP-1-dependent transcription. The transcriptional activity of c-Jun is regulated by phosphorylation at Ser63 and Ser73 through SAPK/JNK (reviewed in 2). Knock-out studies in mice have shown that c-Jun is essential for embryogenesis (3), and subsequent studies have demonstrated roles for c-Jun in various tissues and developmental processes including axon regeneration (4), liver regeneration (5) and T cell development (6). AP-1 regulated genes exert diverse biological functions including cell proliferation, differentiation, and apoptosis, as well as transformation, invasion and metastasis, depending on cell type and context (7-9). Other target genes regulate survival as well as hypoxia and angiogenesis (8,10). c-Jun has emerged as a promising therapeutic target for cancer, vascular remodeling, acute inflammation, as well as rheumatoid arthritis (11,12).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.