

GABARAP Antibody

Catalog No: #32924

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GABARAP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total GABARAP protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human GABARAP.
Target Name	GABARAP
Other Names	MM46; ATG8A; GABARAP-a;
Accession No.	Swiss-Prot:O95166NCBI Gene ID:11337
SDS-PAGE MW	13KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: 1:500 - 1:2000


Immunohistochemistry: 1:50 - 1:200

Immunofluorescence: 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using GABARAP antibody.



Immunofluorescence analysis of MCF-7 cell using GABARAP antibody. Blue: DAPI for nuclear staining.

Background

Gamma-aminobutyric acid A receptors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.