

## Hemopexin Antibody

Catalog No: #35504

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hemopexin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total Hemopexin protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fragment corresponding to a region within amino acids 151 and 438 of Hemopexin.
Target Name	Hemopexin
Other Names	FLJ56652 antibody; HX antibody; HPX antibody; beta-1B-glycoprotein antibody; hemopexin antibody
Accession No.	Swiss-Prot#:P02790;NCBI Gene#:3263
SDS-PAGE MW	52kd
Concentration	0.44mg/ml
Formulation	Rabbit IgG in 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C

## Background

This gene encodes a plasma glycoprotein that binds heme with high affinity. The encoded protein is an acute phase protein that transports heme from the plasma to the liver and may be involved in protecting cells from oxidative stress. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.