

Adiponectin Antibody

Catalog No: #35509

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Adiponectin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total Adiponectin protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 244 of Adiponectin.
Target Name	Adiponectin
Other Names	ACDC antibody; ACRP30 antibody; ADIPQTL1 antibody; ADPN antibody; APM-1 antibody; APM1 antibody; GBP28 antibody; adiponectin antibody; ADIPOQ antibody; adipocyte complement-related 30 kDa protein antibody; 30 kDa adipocyte complement-related protein antib
Accession No.	Swiss-Prot#:Q15848;NCBI Gene#:9370
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Rabbit IgG in 1XPBS, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:3000

Images



Sample (50 ug of whole cell lysate)
A: mouse brown adipose
12% SDS PAGE
#35509 diluted at 1:500

Background

This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The

encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.