

NRTN Antibody

Catalog No: #37199

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NRTN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NRTN protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human neurturin
Target Name	NRTN
Other Names	NTN
Accession No.	Swiss-Prot#: Q99748NCBI Gene ID: 4902Gene Accssion: NP_004549.1
Concentration	2.1mg/ml
Formulation	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human ovarian cancer tissue using #37199 at dilution 1/60.



Immunohistochemical analysis of paraffin-embedded Human brain cancer tissue using #37199 at dilution 1/60.

Background

Neurturin is a member of the TGF-beta subfamily, TRN. This gene signals through RET and a GPI-linked coreceptor, and promotes survival of neuronal populations. A neurturin mutation has been described in a family with Hirschsprung Disease.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.