

ECFP Mouse Polyclonal Antibody

Catalog No: #37991

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ECFP Mouse Polyclonal Antibody
Host Species	Mouse
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IF
Species Reactivity	N/A
Specificity	The ECFP Antibody detects ECFP & CFP tag protein..
Target Name	ECFP
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of ECFP Recombinant protein, using #37991 diluted at 1) 1:5,000 2) 1:10,000

Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.