

Galectin-3 Mouse Monoclonal Antibody

Catalog No: #37995

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Galectin-3 Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6G2
Purification	Affinity purification using immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	Galectin-3 Mouse monoclonal antibody detects endogenous Galectin-3 proteins..
Target Name	Galectin-3
Other Names	Carbohydrate-binding protein 35; CBP 35; CBP35; Gal-3;
Accession No.	Swiss-Prot#:P17931
SDS-PAGE MW	26kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:1000~1:3000

Immunohistochemistry: 1:200

Images



Western blot analysis of Hela cell lysate, using #37995 diluted at 1:3,000.



IHC staining of paraffin-embedded Human thyroid tissue with Galectin-3 mouse mAbB£B"6G2B£B©diluted at 1:200.

Background

Galectin 3 is one of the more extensively studied members of this family and is a 30 kDa protein. Due to a C-terminal carbohydrate binding site, Galectin 3 is capable of binding IgE and mammalian cell surfaces only when homodimerized or homooligomerized. Galectin 3 is normally distributed in epithelia of many organs, in various inflammatory cells, including macrophages, as well as dendritic cells and Kupffer cells. The expression of this lectin is up-regulated during inflammation, cell proliferation, cell differentiation and through trans-activation by viral proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.