

EFHD1 Mouse Monoclonal Antibody

Catalog No: #38043

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	EFHD1 Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB IF
Species Reactivity	Hu Rt Ms
Specificity	EFHD1 Mouse Monoclonal antibody detects endogenous EFHD1 protein.s
Target Name	EFHD1
Other Names	DKFZp781H0842; EF hand domain containing 1; EF hand domain family, member D1; EF-hand domain family, member D1
Accession No.	Swiss-Prot#:Q9BUP0
SDS-PAGE MW	26kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000~1:3000


Immunofluorescence: 1:100-200

Images



Western blot analysis of 1) Mouse spleen tissue, 2) Rat spleen tissue, using #38043 diluted at 1:3,000.

IF analysis of Hela with #38043(Left) and DAPI (Right) diluted at 1:100.



Background

EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation

Note: This product is for in vitro research use only and is not intended for use in humans or animals.