

NUMBL antibody

Catalog No: #38379

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUMBL antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total NUMBL antibody.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human NUMBL.
Target Name	NUMBL
Other Names	CAG3A; CTG3a; NBL; NUMB-R; NUMBLIKE; NUMBR; TNRC23;
Accession No.	Swiss-Prot#: Q9Y6R0NCBI Gene ID: 9253
SDS-PAGE MW	65kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:2000

Images



Western blot analysis of extracts of kidney cell lines, using NUMBL antibody.

Background

Numb-like (Numb or Nbl) binds Notch and is an intracellular modulator of Notch signaling pathway. Numb contains a phosphotyrosine binding (PTB) domain that could bind the Y682 motif present in the cytoplasmic tail of APP. Deduced molecular weight of Numb is 72 kDa.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.