

## DAP antibody

Catalog No: #38736

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	DAP antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB IF
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous level of total DAP antibody.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human DAP.
Target Name	DAP
Accession No.	Swiss-Prot#: P51397NCBI Gene ID: 1611
SDS-PAGE MW	11kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details


Western blotting: □1:500 - 1:2000

Immunofluorescence: □1:20 - 1:100

## Images



Western blot analysis of extracts of various cell lines, using DAP antibody.



Immunofluorescence analysis of MCF-7 cell using DAP antibody. Blue: DAPI for nuclear staining.

## Background

This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.