

β-Tubulin Antibody

Catalog No: #39299



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | β-Tubulin Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Affinity purified |
| Applications | ELISA WB |
| Species Reactivity | Hu Ms Rt |
| Specificity | Not yet tested in other species. |
| Immunogen Type | Peptide-KLH |
| Immunogen Description | nthetic peptide conjugated to KLH |
| Target Name | beta-Tubulin |
| Concentration | 1.0mg/ml |
| Formulation | PBS, pH 7.4 with 0.02% Sodium Azide |
| Storage | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |

Application Details

ELISA WB

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.