

## HA Antibody

Catalog No: #39360

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HA Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Immunogen Type	Peptide-KLH
Immunogen Description	This Polyclonal antibody was raised against a synthetic peptide corresponding to residues (YPYDVPDYA) coupled to KLH.
Target Name	HA
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

ELISA WB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.