

His Antibody

Catalog No: #39361

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	His Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal His-tag antibody was raised against a 6 amino acid peptide. The peptide sequence: HHHHHH
Target Name	His
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.