

ICSBP Antibody

Catalog No: #39367

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ICSBP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal ICSBP antibody was raised against a 14 amino acid peptide at the C terminal residues of human ICSBP protein.
Target Name	ICSBP
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.