

ubiquitin Antibody

Catalog No: #39382

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ubiquitin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	IHC
Species Reactivity	Hu Ms
Specificity	Human,Mus
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide
Target Name	ubiquitin
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded platelet showing cytoplasmic and nuclear staining with anti-ubiquitin antibody at a dilution of 1/100(A). Block with immunizing peptide(B).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.