

Bcl-2 Monoclonal Antibody

Catalog No: #40415



Package Size: #40415-1 50ul #40415-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Bcl-2 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Bcl-2 proteins.
Target Name	Bcl-2
Other Names	BCL2; Apoptosis regulator Bcl-2
Accession No.	Swiss-Prot: P10415NCBI Gene ID: 596
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:1000~2000

Immunohistochemistry: 1:200

Images



Western blot analysis of Hela, diluted at 1:1000



IHC staining of Human tonsil tissue paraffin-embedded, diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.