

## AMACR Monoclonal Antibody

Catalog No: #40421

Package Size: #40421-1 50ul #40421-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	AMACR Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous AMCAR/P504S proteins.
Target Name	AMACR
Other Names	AMACR; Alpha-methylacyl-CoA racemase; 2-methylacyl-CoA racemase
Accession No.	Swiss-Prot: Q9UHK6NCBI Gene ID: 23600
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000

Immunohistochemistry: 1:200

## Images



Western blot analysis of 1) HepG2, 2) Mouse Kidney, diluted at 1:1000.



IHC staining of mouse prostate adenocarcinoma tissue, diluted at 1:200.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.