

## mOrange Monoclonal Antibody

Catalog No: #40480



Package Size: #40480-1 50ul #40480-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	mOrange Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	N/A
Specificity	The antibody detects mOrange & mOrange Tag protein.
Target Name	mOrange
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:5000

## Images



Western blot analysis of recombinant mOrange protein, diluted at 1:10000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.