

PR Monoclonal Antibody

Catalog No: #40481



Package Size: #40481-1 50ul #40481-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PR Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	IHC IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous PR protein.
Immunogen Type	peptide
Immunogen Description	Synthetic Peptide
Target Name	PR
Accession No.	Swiss-Prot: P06401NCBI Gene ID: 5241
SDS-PAGE MW	38,43kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:200 Immunofluorescence: 1:100-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.