

## Survivin Monoclonal Antibody

Catalog No: #40488



Package Size: #40488-1 50ul #40488-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	Survivin Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous Survivin protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	Survivin
Accession No.	Swiss-Prot: O15392NCBI Gene ID: 332
SDS-PAGE MW	66kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-2000

Immunohistochemistry: 1:200-500

## Images



Western blot analysis of 1) HeLa, 2) 293, 3) PC12 using Survivin Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma using Survivin Monoclonal Antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.