

## AChE Polyclonal Antibody

Catalog No: #40541



Package Size: #40541-1 50ul #40541-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AChE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	AChE Polyclonal Antibody detects endogenous levels of AChE protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human AChE.
Target Name	AChE
Other Names	ACHE; Acetylcholinesterase; AChE
Accession No.	Swiss-Prot: P22303NCBI Gene ID: 43
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of MCF7 cells using AChE Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.