

BCoR Polyclonal Antibody

Catalog No: #40644



Package Size: #40644-1 50ul #40644-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BCoR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	BCoR Polyclonal Antibody detects endogenous levels of BCoR protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BCoR.
Target Name	BCoR
Other Names	BCOR; KIAA1575; BCL-6 corepressor; BCoR
Accession No.	Swiss-Prot: Q6W2J9NCBI Gene ID: 54880
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using BCoR Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.