

## CD34 Polyclonal Antibody

Catalog No: #40701



Package Size: #40701-1 50ul #40701-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	CD34 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	CD34 Polyclonal Antibody detects endogenous levels of CD34 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD34.
Target Name	CD34
Other Names	CD34; Hematopoietic progenitor cell antigen CD34; CD antigen CD34
Accession No.	Swiss-Prot: P28906NCBI Gene ID: 947
SDS-PAGE MW	89kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of HepG2 cells using CD34 Polyclonal Antibody



Western Blot analysis of K562 HepG2 cells using CD34 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.