

Choactase Polyclonal Antibody

Catalog No: #40743



Package Size: #40743-1 50ul #40743-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Choactase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Choactase Polyclonal Antibody detects endogenous levels of Choactase protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Choactase.
Target Name	Choactase
Other Names	CHAT; Choline O-acetyltransferase; CHOACTase; ChAT; Choline acetylase
Accession No.	Swiss-Prot: P28329NCBI Gene ID: 1103
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Choactase Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.