

CNPase Polyclonal Antibody

Catalog No: #40768



Package Size: #40768-1 50ul #40768-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CNPase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	CNPase Polyclonal Antibody detects endogenous levels of CNPase protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CNPase.
Target Name	CNPase
Other Names	CNP; 2'; 3'-cyclic-nucleotide 3'-phosphodiesterase; CNP; CNPase
Accession No.	Swiss-Prot: P09543NCBI Gene ID: 1267
SDS-PAGE MW	140kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using CNPase Polyclonal Antibody



Western Blot analysis of various cells using CNPase Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.