

EGF Polyclonal Antibody

Catalog No: #40865



Package Size: #40865-1 50ul #40865-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	EGF Polyclonal Antibody detects endogenous levels of EGF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human EGF.
Target Name	EGF
Other Names	EGF; Pro-epidermal growth factor; EGF
Accession No.	Swiss-Prot: P01133NCBI Gene ID: 1950
SDS-PAGE MW	175kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using EGF Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.