

GFAP Polyclonal Antibody

Catalog No: #40961



Package Size: #40961-1 50ul #40961-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GFAP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	GFAP Polyclonal Antibody detects endogenous levels of GFAP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human GFAP around the non-phosphorylation site of S38.
Target Name	GFAP
Other Names	GFAP; Glial fibrillary acidic protein; GFAP
Accession No.	Swiss-Prot: P14136NCBI Gene ID: 2670
SDS-PAGE MW	95kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of RAT-MUSCLE cells using GFAP Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.