

HSP90 β Polyclonal Antibody

Catalog No: #41041



Package Size: #41041-1 50ul #41041-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HSP90 β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IP IF ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	HSP90 β Polyclonal Antibody detects endogenous levels of HSP90 β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HSP90 β around the non-phosphorylation site of S254.
Target Name	HSP90 β
Other Names	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Accession No.	Swiss-Prot: P08238NCBI Gene ID: 3326
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. Immunoprec


Immunoprecipitation: 2-5 ug/mg lysate.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using HSP90 β Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.