

IRTKS Polyclonal Antibody

Catalog No: #41078



Package Size: #41078-1 50ul #41078-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRTKS Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	IRTKS Polyclonal Antibody detects endogenous levels of IRTKS protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IRTKS.
Target Name	IRTKS
Other Names	BAIAP2L1; IRTKS; Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 1; BAI1-associated protein 2-like protein 1; Insulin receptor tyrosine kinase substrate
Accession No.	Swiss-Prot: Q9UHR4NCBI Gene ID: 55971
SDS-PAGE MW	125kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using IRTKS Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.