

## Musculin Polyclonal Antibody

Catalog No: #41190



Package Size: #41190-1 50ul #41190-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Musculin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Musculin Polyclonal Antibody detects endogenous levels of Musculin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Musculin.
Target Name	Musculin
Other Names	MSC; ABF1; BHLHA22; Musculin; Activated B-cell factor 1; ABF-1; Class A basic helix-loop-helix protein 22; bHLHa22
Accession No.	Swiss-Prot: O60682NCBI Gene ID: 9242
SDS-PAGE MW	165kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of Jurkat cells using Musculin Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.