

NGF Polyclonal Antibody

Catalog No: #41239



Package Size: #41239-1 50ul #41239-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NGF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	NGF Polyclonal Antibody detects endogenous levels of NGF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NGF.
Target Name	NGF
Other Names	NGF; NGFB; Beta-nerve growth factor; Beta-NGF
Accession No.	Swiss-Prot: P01138NCBI Gene ID: 4803
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.


Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using NGF Polyclonal Antibody

Western Blot analysis of various cells using NGF Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.