

RAGE Polyclonal Antibody

Catalog No: #41387



Package Size: #41387-1 50ul #41387-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAGE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	RAGE Polyclonal Antibody detects endogenous levels of RAGE protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human RAGE.
Target Name	RAGE
Other Names	AGER; RAGE; Advanced glycosylation end product-specific receptor; Receptor for advanced glycosylation end products
Accession No.	Swiss-Prot: Q15109NCBI Gene ID: 177
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using RAGE Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.