

## Syk Polyclonal Antibody

Catalog No: #41469



Package Size: #41469-1 50ul #41469-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Syk Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Syk Polyclonal Antibody detects endogenous levels of Syk protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Syk.
Target Name	Syk
Other Names	SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk
Accession No.	Swiss-Prot: P43405NCBI Gene ID: 6850
SDS-PAGE MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of HepG2 cells using Syk Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.