

TWEAK Polyclonal Antibody

Catalog No: #41518

Package Size: #41518-1 50ul #41518-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TWEAK Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	TWEAK Polyclonal Antibody detects endogenous levels of TWEAK protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TWEAK.
Target Name	TWEAK
Other Names	TNFSF12; APO3L; DR3LG; Tumor necrosis factor ligand superfamily member 12; APO3 ligand; TNF-related weak inducer of apoptosis; TWEAK
Accession No.	Swiss-Prot: O43508NCBI Gene ID: 8742
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using TWEAK Polyclonal Antibody



Western Blot analysis of 293T AD293 HELA HepG2 KB cells using TWEAK Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.