

# Tyk 2 Polyclonal Antibody

Catalog No: #41519



Package Size: #41519-1 50ul #41519-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tyk 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Mk
Specificity	Tyk 2 Polyclonal Antibody detects endogenous levels of Tyk 2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Tyk 2 around the non-phosphorylation site of Y1054.
Target Name	Tyk 2
Other Names	TYK2; Non-receptor tyrosine-protein kinase TYK2
Accession No.	Swiss-Prot: P29597NCBI Gene ID: 7297
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293 cells using Tyk 2 Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.