

# XIAP Polyclonal Antibody

Catalog No: #41536



Package Size: #41536-1 50ul #41536-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	XIAP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	XIAP Polyclonal Antibody detects endogenous levels of XIAP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human XIAP around the non-phosphorylation site of S87.
Target Name	XIAP
Other Names	XIAP; API3; BIRC4; IAP3; E3 ubiquitin-protein ligase XIAP; Baculoviral IAP repeat-containing protein 4; IAP-like protein; ILP; hILP; Inhibitor of apoptosis protein 3; IAP-3; hIAP-3; hIAP3; X-linked inhibitor of apoptosis protein; X-linked I
Accession No.	Swiss-Prot: P98170NCBI Gene ID: 331
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.


Not yet tested in other applications.

## Images



Western Blot analysis of 293 cells using XIAP Polyclonal Antibody

Western Blot analysis of 4T1 293 cells using XIAP Polyclonal Antibody



---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.