

GAPDH Polyclonal Antibody

Catalog No: #41549



Package Size: #41549-1 50ul #41549-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GAPDH Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	GAPDH Polyclonal Antibody detects endogenous levels of GAPDH protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GAPDH.
Target Name	GAPDH
Other Names	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH
Accession No.	Swiss-Prot: P04406NCBI Gene ID: 2597
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.04% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 HepG2 MCF7 cells using GAPDH Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.