

## Layilin Polyclonal Antibody

Catalog No: #41630



Package Size: #41630-1 50ul #41630-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Layilin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Layilin Polyclonal Antibody detects endogenous levels of Layilin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Layilin.
Target Name	Layilin
Other Names	LAYN; Layilin
Accession No.	Swiss-Prot: Q6UX15NCBI Gene ID: 143903
SDS-PAGE MW	130kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HuvEc K562 cells using Layilin Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.