

DPYD Polyclonal Antibody

Catalog No: #41635

Package Size: #41635-1 50ul #41635-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DPYD Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	DPYD Polyclonal Antibody detects endogenous levels of DPYD protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human DPYD.
Target Name	DPYD
Other Names	DPYD; Dihydropyrimidine dehydrogenase [NADP(+)]; DHPDHase; DPD; Dihydrothymine dehydrogenase; Dihydrouracil dehydrogenase
Accession No.	Swiss-Prot: Q12882NCBI Gene ID: 1806
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using DPYD Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.