

NTE Polyclonal Antibody

Catalog No: #41751

Package Size: #41751-1 50ul #41751-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NTE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	NTE Polyclonal Antibody detects endogenous levels of NTE protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NTE.
Target Name	NTE
Other Names	PNPLA6; NTE; Neuropathy target esterase; Patatin-like phospholipase domain-containing protein 6
Accession No.	Swiss-Prot: Q8IY17NCBI Gene ID: 10908
SDS-PAGE MW	170kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of AD293 cells using NTE Polyclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.