

## Cryopyrin Polyclonal Antibody

Catalog No: #41768



Package Size: #41768-1 50ul #41768-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cryopyrin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Cryopyrin Polyclonal Antibody detects endogenous levels of Cryopyrin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Cryopyrin.
Target Name	Cryopyrin
Other Names	NLRP3; C1orf7; CIAS1; NALP3; PYPAF1; NACHT, LRR and PYD domains-containing protein 3; Angiotensin/vasopressin receptor AII/AVP-like; Caterpillar protein 1.1
Accession No.	Swiss-Prot: Q96P20NCBI Gene ID: 114548
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using Cryopyrin Polyclonal Antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.