

Ub Polyclonal Antibody

Catalog No: #41812



Package Size: #41812-1 50ul #41812-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ub Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Ub Polyclonal Antibody detects endogenous levels of Ub protein only when non-acetylation at K63.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Ub around the non-acetylation site of K63.
Target Name	Ub
Other Names	UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2
Accession No.	Swiss-Prot: P62987/P62979/P0CG47/P0CG48NCBI Gene ID:
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Ub Polyclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.