

PTN Polyclonal Antibody

Catalog No: #41889



Package Size: #41889-1 50ul #41889-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	PTN Polyclonal Antibody detects endogenous levels of PTN protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PTN.
Target Name	PTN
Other Names	PTN; HBNF1; NEGF1; Pleiotrophin; PTN; Heparin-binding brain mitogen; HBBM; Heparin-binding growth factor 8; HBGF-8; Heparin-binding growth-associated molecule; HB-GAM; Heparin-binding neurite outgrowth-promoting factor 1; HBNF-1; Osteoblast-specific facto
Accession No.	Swiss-Prot: P21246NCBI Gene ID: 5764
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using PTN Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.