

TNAP Polyclonal Antibody

Catalog No: #41929



Package Size: #41929-1 50ul #41929-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TNAP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Rt
Specificity	TNAP Polyclonal Antibody detects endogenous levels of TNAP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human TNAP.
Target Name	TNAP
Other Names	ALPL; Alkaline phosphatase, tissue-nonspecific isozyme; AP-TNAP; TNSALP; Alkaline phosphatase liver/bone/kidney isozyme
Accession No.	Swiss-Prot: P05186NCBI Gene ID: 249
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using TNAP Polyclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.