

Probetacellulin Polyclonal Antibody

Catalog No: #42006



Package Size: #42006-1 50ul #42006-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Probetacellulin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Probetacellulin Polyclonal Antibody detects endogenous levels of Probetacellulin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Probetacellulin.
Target Name	Probetacellulin
Other Names	BTC; Probetacellulin
Accession No.	Swiss-Prot: P35070NCBI Gene ID: 685
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of KB cells using Probetacellulin Polyclonal Antibody

Western Blot analysis of KB cells using Probetacellulin Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.