

GFP-Tag Rabbit Polyclonal Antibody

Catalog No: #T508



Package Size: #T508-1 50ul #T508-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GFP-Tag Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Tag
Specificity	GFP-Tag Rabbit Polyclonal antibody detects endogenous levels of C-terminal, internal, and N-terminal GFP-tagged proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant GFP expressed in E.coli
Target Name	GFP-Tag
SDS-PAGE MW	N/A
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Application Details

WB: 1:3000~1:10000

Images



Western blot analysis of PET-28a GFP recombinant protein expressed in E.coli, using GFP-Tag Rabbit Polyclonal antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.