

GST-Tag Mouse Monoclonal Antibody

Catalog No: #T509



Package Size: #T509-1 50ul #T509-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GST-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Tag
Specificity	GST-Tag Mouse Monoclonal antibody detects endogenous levels of C-terminal, internal, and N-terminal GST-tagged proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Full-length GST protein
Target Name	GST-Tag
SDS-PAGE MW	N/A
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C. Stable for 12 months from date of receipt

Application Details

WB: 1:3000~1:10000

Images



Western blot analysis of GST-tag fusion protein, using GST-Tag Mouse Monoclonal antibody (lane 1: 1:2000; lane 2: 1:5000).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.