

RFP Mouse Polyclonal Antibody

Catalog No: #T511

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RFP Mouse Polyclonal Antibody
Host Species	Mouse
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IP
Species Reactivity	N/A
Specificity	The RFP tag antibody can recognize RFP tag fusion proteins.
Target Name	RFP
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of RFP Recombinant protein, using #T511 diluted at 1) 1:5,000, 2) 1:10,000.

Background

The red fluorescent protein cloned from *Discosoma coral* (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea jellyfish*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.